


ALL INFORMATION CONTAINED IN THIS DOCUMENT SHOULD BE USED FOR REFERENCE PURPOSES ONLY & ARE SUBJECT TO CHANGE WITHOUT PRIOR NOTIFICATION. ALSO, THIS DOCUMENT IS CONSIDERED PROPRIETARY & THE SOLE PROPERTY OF CRYOMECH, INC. & SHOULD NOT BE REPRODUCED IN ANY MANNER WITHOUT PRIOR WRITTEN APPROVAL.

B

B



CP2870 COMPRESSOR PACKAGE (WATER COOLED)

POWER CORD

COLD HEAD MOTOR CORD (SAME LENGTH AS FLEX LINES)

(2) HELIUM FLEX LINES (1) 1/2" ID (1) 3/4" ID  
STANDARD LENGTH: 10ft  
LONGER LENGTHS AVAILABLE


3/4" AEROQUIP

3/4" AEROQUIP

1/2" AEROQUIP

3/4" AEROQUIP

NITROGEN FLEX LINE 1/4" ID  
STANDARD LENGTH: 10ft  
LONGER LENGTHS AVAILABLE




NITROGEN GENERATOR ASSEMBLY (STAND OPTIONAL)



210 LITER LIQUID NITROGEN DEWAR ASSEMBLY

AL300 COLD HEAD



AIR COMPRESSOR (OPTIONAL)

A


A



CP2870 COMPRESSOR PACKAGE (AIR COOLED)

POWER CORD

COLD HEAD MOTOR CORD (SAME LENGTH AS FLEX LINES)



MANUAL




TOOL KIT



1 YEAR SPARE PARTS

ECN	ZONE	REV	DESCRIPTION	DATE
3072	-	C	ADDED CP2800 COMPRESSOR PACKAGE	15MAY13
2682	-	B	UPDATED COMPRESSOR SIDE PANELS	23FEB12
1419	-	A	SEE ECN 1419	03JUL07

THIRD ANGLE PROJECTION



**CRYOMECH, INC**

113 FALSO DRIVE SYRACUSE, N.Y. 13211  
Tel: (315)455-2555 Fax:(315)455-2544

ALL DIMENSIONS ARE IN INCHES  
TOLERANCES: FRACTIONS ±1/32 DECIMALS XXX=±.005 XX=±.01 X=±.1  
ANGLES ±.5° ALLOVER FINISH 25 Ra  
CONCENTRICITY .005 TOTAL INDICATOR RUN OUT  
PERPENDICULAR ±.002 (UNLESS OTHERWISE SPECIFIED)

NAME: LNP60 LIQUID NITROGEN PLANT

DWG #: LNP60 SYSTEM SHT 1 OF 1

MATERIAL: AS STATED

DWN BY: BL DATE: 18AUG04 SCALE: TO FIT DWG SIZE: A

1ST CHK: RED DATE: 26AUG04 2ND CHK: NLO DATE: 26AUG04

CUST: BTA: